

---

# MAKING THE IPAD ACCESSIBLE...WHO KNEW?

---

Oklahoma Abletech AT Expo

April 4, 2014

Tulsa, OK

Dr. Sandra Wright Ph.D., CCC-SLP, BCS-CL

---

---

# AN INVALUABLE THERAPEUTIC TOOL

---

- Many communication applications
    - proloquo2go
    - soundingboard
    - go talk app
    - sonoflex
    - verbally
    - many many more
-

---

# DOWNFALLS OF USING AN IPAD

---

- durability
  - multitasking
  - so many apps so little time
  - accessibility
-

---

Times have changed though!!! At least for  
access

THE IPAD HAS HEAD  
CONTROL ACCESS  
BUILT-IN AND SCREEN  
AS A SWITCH ACCESS

---

---

# HOW TO SET UP SWITCH CONTROL

---

- It is built into ios7 (a perk for the update)
  - go to settings
  - go to accessibility
  - go to switch control
-

---

# CAVEATS

---

- some apps don't have switch control "built in"
  - some apps have switch control setting options but they can only be managed within the app itself not from the ios setting button
  - if using your screen as a switch, no escape out of this, so be certain this is what you want...
-

---

# SETTING UP EXTERNAL SWITCHES

---

- Direct connect a switch (using the same port as you would to charge your device or sync it to a computer)
  - options of (one switch control up to 5 switch control)
  - Bluetooth connect a switch (several options available)
  
  - Most switches are "plug and play"
-