

Make 'n Take Sheet

Discover How Easy It Is To Make Your Own Assistive Technology!

Positioning Tool

For the children that need their legs propped up a little bit to make seats more comfortable.

Materials:

Bungee Cord
Pool Floatie
Scissors



Make-N-Take Price: \$1.00-\$10.00
(depending on if you already had supplies)

Directions:

1. Cut the pool floatie to desired length (about the width of a carseat or stroller).
2. Thread bungee cord through pool floatie.
3. Attach to carseat or stroller.

Helpful Hints: You can get a thinner floatie or a thicker floatie, depending on the amount of support needed.

Oklahoma ABLE Tech, the statewide Assistive Technology Act Program, is funded through the Administration for Community Living of the U.S. Department of Health and Human Services, and maintains coordination and collaboration efforts with partners throughout the State of Oklahoma.

To contact Oklahoma ABLE Tech call 1-800-257-1705 or INFO-line 1-888-885-5588