

Utensil Grip

The modified utensil grip enables people with low muscle tone or other physical disabilities to grasp a utensil. This can also be useful with toothbrushes, combs, pencils, paintbrushes and etc.

Materials:

- 1 – Utensil
- 1 – Hair Roller or any type of sponge tubing (approx. 1 inch diameter)

Directions:

1. Remove plastic part of curler
2. Place Styrofoam curler or foam tubing over the handle of the utensil



Cost Comparison:

Market Price: \$2.00

Make-N-Take Price: \$.50